
Project Folders and Files Structure:



Used Packages:

- 1. "react-router-dom": "react-router" package provides the core routing functionality for React Router. Instead of installing it directly, I have installed "react-router-dom" and It installs "react-router" as a dependency. In my web app there are two types of routes, public routes and private routes and is declared in routes.tsx file.
- 2. "axios": "axios" is a library to make XMLHttpRequests from the browser. "axios" has wide browser support, client side support for protecting against XSRF, it automatically transforms to JSON data, it supports the promise API.
- 3. "universal-cookie": "universal-cookie" is the most popular package for managing cookies. [1]




fig 1: Downloads statistics

And it's easy to use and well documented.

- 4. "react-select-async-paginate": Used to create TypeId select field, that supports pagination on options menu scroll.
- 5. Material UI ("@material-ui/core", "@material-ui/icons", "@material-ui/lab"): Material-UI is a popular React UI framework and an open-source project that features React components which implement Google's Material Design.

UI Design and Functionality:

Login Page (Desktop View): Without logged in, user will be redirected to login page for any routes.




Fig2: Login page (Desktop View)

Login Page (Mobile View: Model Nexus 4 - Google Chrome emulated device)



Fig3: Login Page (phone view)

Available Medical Devices Page (Desktop View): In this page user can select how many device model will be shown. And this page has a pagination. If user clicks on ADD NEW DEVICE MODEL button then he will see a modal dialog with form to create New Device Modal.




fig4: Available Medical Devices Page (Desktop View)




fig5: Available Medical Devices Pagination (Desktop View)

Available Medical Devices Page (Phone View):






fig6: Available Medical Devices (Phone view)

fig7: Available Medical Devices Pagination (Phone View)

Model Data Overlay Dialog (Desktop View): After click on a model device, user will see the model device data in an overlay Dialog. Before data is fully loaded user will see a circle progress bar on the dialog.




fig8: Model Data Dialog (If any Data found)





fig9: Model Data Dialog (If no data found)

Model Data Overlay Dialog (Phone View):






fig10-11: Model Data Dialog- Phone View (left- If data found, right-If no data found)

Add New Device Model (Desktop View): Create error message or success message will be shown above first input field.




fig12: Add new device model form (Desktop View)

Add New Device Model (Phone View):



Fig13: Add new device model form with dialog (Phone view)

Thank you